

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Tetra-n-butylammonium hydroxide - (20% solution in water) for synthesis

Product Number:

8.18759

TEST	Specification
------	---------------

Appearance (Color)

according to color reference solution Ph.Eur.:

not more intensely colored than reference standard Y5

Complying

Assay

(acidimetric)

18.0 - 22.0 %

Remarks:

Due to its specific melting range the product may be solid, liquid, a solidified melt or a supercooled melt.

Specification: PRE.0.ZQ5.10000025296

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

