

HPLC PEEK NUTS / FERRULES



SPECIFICATION

Accuracy: Ensures leak-free high-precision connections

Temperature Range: -20C to 120C (continuous use)

Model No: Custom/Universal Fit for HPLC

Core Components: PEEK body with precision-machined threading

Feature: Chemically inert, high pressure resistance, reusable, biocompatible

Equipment Materials: PEEK (Polyether Ether Ketone)

Usage: Used to connect tubing in HPLC systems

Dimension (L*W*H): 1/4-28 thread (Standard), 1.6 mm to 1/16 inch OD tubing compatible