

Distillation Unit SS 5L/hr



Horizontal thick wall heat resistant borosilicate glass boiler fitted with 3 KW chromium plated heater.

High efficiency coil condenser ensure low temperature, high output quality distillate of conductivity 3 to 4 μ s/cm.

Compact construction suitable for bench and wall mounting.

Screwthread connectors enable easy fitting or removal of hose safely.

Built-in constant level device with Vendura-Nova valve for drainage provides easy descaling without dismantling.

Safety thermostat built into the heater protect the water still in the event of water supply failure and automatically resets when the water supply is restored

Distillate output	4 litres per hour, single distilled
pH	5.0 to 6.5
Conductivity	3.0 to 4.0 μ s/cm
Temperature	35 to 50 $^{\circ}$ C
Pyrogen Content	Pyrogen free
Power Requirement	220/240v , 50/60 Hz single phase, 13A fused
Dimension (d x w x h)	50 x 15 x45 cm
Catalogue Number	P21001